

Slides for Pediatric Model

Preschool Children
The “E” points to a Bird, a Rabbit, a Girl or a Boy.




Test No. 1, SMPS-B2 – Binocularity
Tests performance of both eyes together.




**Test No. 2, SMPS-30-1 – Acuity
Right and Left Eye, 20/30**
The “E” points in a different direction for each eye.




**Test No. 3, SMPS-30-2 – Acuity
Right and Left Eye, 20/30**
The “E” points in a different direction for each eye.




**Test No. 4, SMPS-30-3 – Acuity
Right and Left Eye, 20/30**
For each eye, all the “E”’s point to a different object.
Also available 20/40 acuity.



**Test No. 5, SAS-2 – Acuity
Right and Left Eye**
Conventional Tumbling “E” symbol. Child states which way the “E” points in each square. Acuity levels of 20/20, 20/30, 20/40 and 20/50.



Test No. 6, SP3-3 – Heterophorias
Muscle imbalance – vertical and lateral – is tested at predetermined pass/fail levels. This test is administered first at far, then at near.



Test No. 7, SPCD-1 – Color Perception
Titmus Pseudo-Isochromatic Plates reliably detect gross color anomalies, but do not classify the type of deficiency.



Test No. 8, SA-BRL-4 – Visual Acuity
Results calibrated at 20/20, 20/30, 20/40, 20/50, 20/70, 20/100 and 20/200. Tests can be administered binocularly or with either eye occluded.

Improvements have been made to test slides to enhance viewing performance. Patent Pending.

Toll Free 800-446-1802